

# Pacemaker Leads

## Sensing Performing Factors

### Factors

- Sense amplifier
- P and R wave amplitude
- Tissue - electrode proximity
- Depolarization wave factor
- Electrode surface area

### Influenced by:

- IPG
- Patient
- Lead design & placement
- Lead placement
- Electrode design

# Better Sensing with Selox ST/JT

## Tissue - Electrode Proximity

Optimal myocardial contact is crucial for signal amplitude detection and slew rate

- Selox has an optimized lead design of tip & tines  
**Benefit: Optimal myocardial contact** (atrium and ventricle)
- Short bipole in Selox leads (15 mm)  
**Benefit: Assures excellent signal specificity and high slew rates**

# Geometry of the Tines

## Tissue - Electrode Proximity

Selox has 3 tines

Maximum width of each tine is 1.6 mm for Selox

### Selox ST/JT

Length

4 mm


Geometry

Structure more  
shapely

# Geometry of the Tines



## More Reliable Fixation with Selox



# Sensing with Selox ST/JT

## Electrode Surface Area

Ideal leads are compromised of small geometric and large electrically active surfaces.

Selox feature fractal coating<sup>1,2,3</sup> of BIOTRONIK

### Benefits:


- Sensing of smallest amplitudes
- Superb signal to noise ratio
- Optimal compatibility to automatic capture control algorithms<sup>3</sup>

<sup>1</sup> Israel et al., Herzschrittmacher 18(1), 1998

<sup>2</sup> Bolz et al., Herzschrittmachertherapie und Electrophysiologie 7 (4)

<sup>3</sup> Philos DR Active Capture Control, IDE No. G010286 Clinical Report (2002)

# Fractal Coating



The "Cauliflower" structure  
of a fractal coated electrode tip  
(self repeating structure)

Small geometric size: efficient  
transfer of energy from tip  
electrode to heart tissue

Large bio-effective surface area for  
the clear detection of intrinsic  
cardiac signals by the pacemaker  
and elimination of polarization  
artifacts

# Pacemaker Leads

## Pacing Performing Factors

### Factors

- Myocardial tissue at pacing position

- Small electrode size

- Acute inflammatory response

- Scar tissue due to enduring inflammatory processes

### Influenced by:

- Patient

- Lead design

- Lead design

- Lead design

# High Pacing Impedance with Selox

## Small Electrode Size with fractal coating

- 1.3 mm<sup>2</sup> electrode tip (Selox)

### Benefits:

- High electrical density and low thresholds
- High impedance



# Optimal Pacing with Selox

## Reduction of acute inflammatory responses via

Selox has a Platinum-Iridium lead tip

### Benefit:

- Excellent biocompatibility for minimum inflammatory responses

Selox has fractal coating

### Benefit:

- Optimal charge transfer between electrode and myocardium
- Long-term stability of pacing threshold

# Lower Acute Pacing Thresholds with Selox ST/JT

## Reduction of acute inflammatory responses via

Selox leads have 0.75 mg Steroid near the lead tip


### Benefit:

- Prevents tissue inflammation
- No post-implant threshold raises

# Selox ST/JT – Clinical Evaluation:




## Sensing [mV] Atrium



# Selox ST/JT – Clinical Evaluation:



## Sensing [mV] Ventricle




Impl: implant; PHD: pre-hospital discharge; 1-M: 1 month after discharge, 3-M: 3 months

# Selox ST/JT – Clinical Evaluation:



## Threshold [V] Atrium



# Selox ST – Clinical Evaluation:



## Threshold [V] Ventricle

